

SL NO	TEST CODE	INVESTIGATION NAME	NOTE	TEMPERATURE	REPORTING TIME	COST
<b>SEROLOGY</b>						
1.	HIVR	HIV RAPID, CARD METHOD	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
2.	HIV	HIV 1 & 2 ANTIBODIES, SCREENING TEST	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM. PATIENT CONSENT ON HIV CONSENT FORM IS MANDATORY.	2-8 <sup>0</sup> C	SAME DAY	
3.	HBSAG	HEPATITIS B SURFACE ANTIGEN (HBsAg), QUANTITATIVE	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
4.			For other advanced tests referred out pls refer serial number 179-183			
5.	HCV	HEPATITIS C ANTIBODY (ANTI-HCV)	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
6.	WIDAL	WIDAL SLIDE AGGLUTINATION TEST	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
7.	WIDALT	WIDAL TUBE AGGLUTINATION TEST	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
8.	VDRL RPR	VDRL/RPR SLIDE FLOCCULATION	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
9.	DENGUE- Ab	DENGUE FEVER ANTIBODIES PANEL, IgG & IgM	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
10.	DENGUEP NS1, Ab	DENGUE FEVER COMBINED PANEL (DENGUE FEVER ANTIBODIES IgG & IgM AND DENGUE FEVER NS1 ANTIGEN)	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
11.	DENGUE-Ag	DENGUE FEVER NS1 ANTIGEN, RAPID	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
12.	HPAG	H-PYLORI ANTIGEN	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	
13.	VDRL	VDRL CARD METHOD	COLLECT 2 ML BLOOD IN A YELLOW TOP (SST) TUBE OR 1 ML SERUM.	2-8 <sup>0</sup> C	SAME DAY	