

SL NO	TEST CODE	INVESTIGATION NAME	NOTE	TEMPERATURE	REPORTING TIME	COST
CLINICAL PATHOLOGY						
1.	URINER	URINE EXAMINATION, ROUTINE/MICROSCOPY	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	
2.	RUS	RANDOM URINE SUGAR	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	
3.	FUS	FUS: FASTING URINE SUGAR	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER. OVERNIGHT FASTING IS MANDATORY	2-8 ⁰ C	SAME DAY	
4.	PPUS	PPUS: POST PRANDIAL URINE SUGAR	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER. COLLECT PPUS SAMPLE EXACTLY 2 HOURS AFTER STARTING MEAL.	2-8 ⁰ C	SAME DAY	
5.	BILE	BILE SALT, URINE	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	
6.	BILEP	BILE PIGMENTS, URINE	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	
7.	UACETONE	KETONE BODIES, URINE	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	
8.	UALB	ALBUMIN, URINE	COLLECT 10 ML (5 ML MIN) IN A SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	
9.	STOOLR	STOOL EXAMINATION, ROUTINE/MICROSCOPY	5g (3g MIN) STOOL IN A LEAK PROOF SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	

10.	STOOLOB	STOOL EXAMINATION, OCCULT BLOOD	5g (3g MIN) STOOL IN A LEAK PROOF SCREW CAPPED CONTAINER. PATIENT SHOULD AVOID MEAT PRODUCTS, PAIN KILLERS, IRON TABLETS, VITAMIN C (ASCORBIC ACID) AND HARD TOOTH BRUSH FOR 3 DAYS PRIOR TO SPECIMEN COLLECTION.	2-8 ⁰ C	SAME DAY	
11.	STOOLRS	STOOL EXAMINATION, REDUCING SUBSTANCES	5g (3g MIN) STOOL IN A LEAK PROOF SCREW CAPPED CONTAINER.	2-8 ⁰ C	SAME DAY	
12.	SHD	STOOL EXAMINATION, HANGING DROP PREPARATION	5g (3g MIN) STOOL IN A LEAK PROOF SCREW CAPPED CONTAINER. SPECIMEN MUST REACH LAB WITHIN 1 HOUR OF COLLECTION AT ROOM TEMPERATURE. Done for Cholera	18-22 ⁰ C	SAME DAY	
13.	SEMEN	SEMEN ANALYSIS	SUBMIT SEMEN IN A SCREW CAPPED STERILE CONTAINER. 3 DAYS ABSTINENCE REQUIRED SAMPLE SHOULD REACH THE LAB WITHIN ONE HOUR	18-22 ⁰ C	SAME DAY	